Pall Corporation

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1. INTRODUCTION

The following procedures must be followed for the installation of SUPRADisc modules.

These instructions as well as the product data sheet should be read thoroughly as they contain valuable information gained by extensive experience. It is very important that all instructions are carefully followed and, where appropriate, should be incorporated into the operator's standard operating procedures.

If some of the procedures do not suit your needs, please consult Pall or your local distributor before finalizing your system.

Use of this product in a manner other than in accordance with Pall's current recommendations may lead to injury or loss. Pall cannot accept liability for such injury or loss.

2. SPECIFICATIONS

The maximum working pressure and temperature may vary between filter style and filter media. Please check the appropriate data sheet for details or contact Pall or your local distributor.

Operation outside the specifications and with fluids incompatible with construction materials may cause personal injury and result in damage to the equipment. Incompatible fluids are those which chemically attack, soften, stress attack or adversely affect the materials of construction.

Please refer to the Chemical Compatibility Guide USTR 2305.

3. RECEIPT OF EQUIPMENT

- Store the filter in clean, well vented, odorless and dry conditions between 0°C (32°F) and 45°C (110°F) without exposure to irradiation sources like direct sunlight, and wherever practical in the packaging as supplied.
- DO NOT remove from packaging until prior to installation.
- Check that the bag and packaging are undamaged prior to use.
- Ensure that the type of filter selected is suitable for the application. In addition to the part number, each filter assembly is identified by a unique identification batch number.
- The SUPRADisc modules should be used within a period of 36 months after production. Exception: BIO SD modules must already be used within 18 months after production.

4. INSTALLATION AND OPERATION

SUPRADisc modules are high quality products manufactured to exacting standards. It is essential to take care when handling and installing them into filter housings.

Before installation, it is essential to verify that the filter type selected is suitable for the fluid to be filtered, and to follow the instructions listed below.

- It is recommended to wear gloves to prevent accidental contamination of the filter.
- The installation and operating instructions supplied with the filter housing should be followed carefully.
- For a safe installation of the modules it is recommended to use the corresponding lifting/inserting aid.
- The housing must be completely de-aerated
- Open the filter units only in a pressureless condition.
- SUPRADisc Modules contain nominally rated depth filter sheet media. For processes where microbiological stability is the target, it is recommended to use membrane cartridge filters qualified for microorganism reduction or pasteurization as final process step.

4.1 Double O-ring style SUPRADisc modules

Ensure that O-rings are undamaged and are correctly positioned in the groove.

Check that the sealing surface on the filter housing is clean and undamaged.

To assist ease of fitting, it is strongly recommended that O-rings be lubricated by dipping the gaskets in a suitable liquid which is compatible with the fluid to be filtered.

Water with the same quality as that used for final rinsing of the installation is a satisfactory lubricant in many cases. For advice on other lubricants, please contact Pall.

Low boiling point lubricants (e.g. ethyl or isopropyl alcohol) must not be used if the installed filter is to be subsequently steam sterilized or exposed to temperatures above the boiling point of the lubricant. The high vapor pressures between the O-rings under these conditions can result in damage to the O-ring adaptor.

Insert the SUPRADisc modules with a gentle twisting motion to assist wetting of the surfaces. Gently ease into place. Do not attempt to force the modules into position.

Finally, twist the modules clockwise to engage the retaining lugs within the filter head.

4.2 Flat gasket style SUPRADisc modules

Ensure that gaskets are undamaged and correctly fitted into modules grooves at each end.

Check that sealing faces on housing and seal nuts are clean and undamaged.

No wetting of gaskets is required.

Slide SUPRADisc modules over tie rod and secure in place with seal nut and HAND TIGHTEN. Then turn the tightening sleeve 3 times per installed module.

5. RINSING

Depending on the application, rinsing with cold or warm water in a forward flow direction is recommended prior to putting the SUPRADisc modules into operation. In special cases rinsing can be carried out with product or a product-compatible liquid. The recommended rinsing volume is 50 L/m² of filter area (for SUPRADisc HP modules 100L/m²).

The recommended flow rate for the rinsing cycle is 1.5 times the process flow rate with a maximum temperature of 85°C (185°F).

6. REGENERATION

SUPRADisc 1:

SUPRADisc 1 Modules can only be used in the forward flow direction. This applies to all processes including rinsing, regeneration, sanitization and filtration. Reverse mechanical impact causes damage to the filter sheet media. Pressure spikes from downstream of the filter or Vacuum after steaming must be avoided. The use of check valves downstream (on the filtrate side) of the housing is recommended.

SUPRADisc II:

SUPRADisc II Modules are designed to withstand flow in both directions. Regeneration of the sheet media may be enhanced with a backwash procedure. If reverse flush procedures are employed, Pall recommends the use of a backwash device. Please contact Pall for ordering information for the backwash device.

When the backwash devices are installed in combination with the SUPRADisc II modules, it is possible to cold water rinse the modules in the reverse direction up to a maximum of 0.5 bar dp. To avoid pressure shocks, check valves are recommended downstream (on the filtrate side) of the housing.

6.1 Steam in place

Please refer to the appropriate Pall product information literature for products which can be steam sterilized in place. Detailed sterilization

procedures can be found in Pall publication USD 2299.

Product and back pressure surges must be avoided during flushing, steaming or filtration in order to avoid damage and to achieve optimal use of the SUPRADisc modules.

In order to avoid injuries from escaping steam or hot water, it is imperative to observe the corresponding safety and protection measures.

7. EC DIRECTIVE 94/9/EC (EQUIPMENT FOR USE IN POTENTIALLY EXPLOSIVE ATMOSPHERES, ATEX)

SUPRADisc modules comply with the ATEX directive when installed in a compliant Pall filter housing or assembly, but are not in themselves required to be labeled with the ATEX marking. Under the terms of the directive, filter modules are not considered to be equipment capable of autonomous function, but may be thought of as components being essential for the operation of the equipment. As such, the conformity of the filter modules has been assessed as an integral part of the overall assembly.

7.2 Hot water Sanitization

SUPRADisc Modules can be sanitized with hot water. This sanitization is best done following the hot water rinsing step described in 5.

The hot water temperature should be adjusted to > 85°C on the outlet pipe of the housing. When 85°C is reached, circulate hot water for 25 - 30 minutes.

8. FILTER MODULE REPLACEMENT

SUPRADisc modules should be replaced according to the GMP requirements of the process. Where filters are used for more than one manufacturing batch, replacements are recommended when the maximum allowable differential pressure has been reached, (refer to appropriate Pall datasheet) or when the flow rate has become unacceptable, whichever occurs first. Discard filter in accordance with local Health and Safety and Environmental Procedures. No attempt should be made to clean disposable filter modules.

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STEAM STERILIZATION OF PALL® SUPRADISC™ MODULES (USD2299)

Important Notice

Refer to safety instructions before use. Safety instructions in this language are available from Pall.

1. IMPORTANT RECOMMENDATIONS

This section includes important procedures and guidelines; it should be reviewed carefully before implementing steam sterilization protocols.

This publication describes important aspects of the procedures to be adopted when steam sterilizing Pall filter assemblies. It cannot take account of particular features of individual systems. Should you find difficulty in applying these recommendations or have any questions concerning steam sterilization in general, please contact your nearest Pall office.

Please note that Pall SUPRADisc filters are not supplied sterile and validation of any sterilization procedure is the responsibility of the user.

1.1 Installation

The filter assembly should be installed in such a manner that condensate from the steam supply cannot accumulate in the housing and that the open end of the assembled filter module is orientated downwards. It is preferable that the pipework downstream of the filter assembly is kept as short as possible. Critical pipe lengths are shown in the procedure diagrams.

1.2 Control of Steam Sterilization

Steam used for sterilization must be saturated and free from condensation. Superheated steam must not be used. Introduction of steam into the system should be in such a way as to prevent 'air traps' forming. Air pockets can inhibit steam flow and produce regions where inadequate sterilizing conditions are achieved. This requires special attention where steam is introduced from more than one position. Adequate means condensate drainage should be employed to ensure that steam is free from condensate.

Condensate will wet hydrophilic filter assemblies, increase differential and reduce steam flow.

It is important to consider:

- Steam supply
- System to be steamed (adequate drains etc)
- Pipe orientation
- Pipe insulation

Steam and air pressure should be regulated carefully to avoid over-pressurization and damage to modules. Accurate and calibrated pressure gauges are important. Pressure differential should be kept to a minimum, but it should not

exceed 300 mbar (4.3 psi) in flow direction. At the completion of steam sterilization air should be introduced to replace the steam; compensation for steam collapse is important to prevent a vacuum forming, which may cause filter damage, leakage via pressure seals, or vessel collapse. Pressure in the reverse flow direction must be avoided.

1.3 Direction of Steam Flow

Depth filter modules are intended to be steam sterilized in situ by steam flow in the normal flow direction (out to in)

1.5 Steam for Sterilization

The steam should be free of particulate matter such as rust and pipe scale as these will be removed by the filter to be sterilized and shorten its life. Pall porous stainless steel filters are suitable for the filtration of steam and appropriate assemblies may be obtained from Pall.

1.5 Flushing

It is strongly recommended that the process filter assembly and associated downstream equipment is flushed after steam sterilization and before filtration to remove any residue originating from the steam and trace amounts of filter extractables remaining after sterilization. Flushing of the filter modules should be carried out in the normal direction with a minimum flush volume of 50 L/m² filter area is recommended. (For SUPRADisc HP modules 100 L/m²) The flow rate should be equivalent to 1.5 times the process flow rate.

2. IN SITU STEAM STERILIZATION PROCEDURES

Pressure Gauges

Pressure gauges which can be read with accuracy over the range of 0 - 3 bar (0 - 43.5 psi) must be installed to monitor steam pressure and differential pressure across the filter assembly during the sterilization cycle.

Steam Sterilization Conditions (see also figure 1)

To ensure effective sterilization, steam temperature (measured at position T) in the assembly should be held at a minimum of 121 °C (250 °F) (1.1 barg (15.9 psi) saturated steam) for the minimum time validated by the user as necessary to achieve system sterilization. The maximum temperature must not exceed 125 °C (257 °F) at a pressure of 1.32 bar (19.1 psi)

Procedure

- Ensure that all valves are closed.
- Fully open drain valve C.
- Fully open condensate drain valve I, open condensate drain valve J, open venting valve G.

- Adjust steam sterilization pressure P4 to 300 mbar (4.3 psi) above the saturated steam pressure P2 necessary for sterilization of the filter. After condensate has been expelled from I, partially close valve I.
- Slowly open steam valve B to admit steam to system. After condensate has been expelled from J, partially close valve J.
- Partially close vent valve G when steam flow is evident, ensuring that pressure at P2 remains within 300 mbar (4.3 psi) of pressure at P1. Partially open drain valve D to drain condensate and partially open valve E.
- Permit steam to flow through the system until steam pressure is stabilized, and adjust the regulated steam supply until the validated temperature is achieved at position T. Monitor temperature at T for the necessary sterilization time. Ensure that pressure at P2 remains within 300 mbar (4.3 psi) of pressure at P1.
- When sterilization time is complete, close valve B. Valve D remains partially open, close valve I completely. When P2 is pressureless, then open valve G and drain valve J completely. Close valves D and E. When the gas pressure has dropped at position T (temperature = 100 °C [212 °F]), open valve A and flush the system with 150 L/m² of ambient temperature water for 10 minutes. After flushing close valve G and drain valve J partially and completely open valve D. Filter assembly is now ready for use.

3. GENERAL GUIDELINES FOR STEAM STERILIZATION OF FILTER ASSEMBLIES

Important:

It is the responsibility of the user to validate the effectiveness and safety of procedures used to steam sterilize process equipment and filter assemblies. The following guidelines are intended only to highlight some aspects of such procedures which require special attention. For further assistance or information, please contact Pall Scientific and Laboratory Services.

3.1 Filter Sizing and Steam Supply

Filter modules should be sized appropriately for product filtration, any gas or air flow and to permit adequate steam flow to sterilize effectively the downstream equipment. Failure to take account of steam flow requirements may result in filter damage, caused by high differential pressures at elevated temperatures, and possible non-sterility of downstream equipment.

3.2 Differential Pressure

During steam sterilization of downstream equipment, differential pressure across the filter assemblies must not exceed 300 mbar (4.3 psi) in the forward direction. Reverse pres-

surization must be avoided in order to prevent filter damage.

3.3 Monitoring of Temperature and pressure

It is important to monitor temperature and pressure in downstream equipment to ensure that:

- Validated sterilizing conditions have been achieved
- Excessive differential pressures are not experienced across filter assemblies
- A sudden fall in pressure due to steam collapse does not compromise downstream equipment

Caution:

Where vessels in the downstream systems are unable to withstand negative pressure without collapse, appropriate safety devices must be fitted.

Product and pressure surges must be avoided to prevent damage and to allow optimal use of the Pall SUPRADisc modules.

3.4 Air Entrapment

It is important to ensure that valve sequences do not lead to entrapment of air pockets in the process equipment, as sterility may be compromised.

3.5 Condensate Drainage

Adequate means for condensate drainage should be employed to ensure that steam is free from condensate. Condensate in the modules can wet the depth filter media, increase differential pressures across the filters and reduce steam flow. Provision should be made for drainage of condensate from process equipment following steam sterilization where such condensate is undesirable for operational reasons.


4. Safety Instructions

- Steaming leads to a heating of the Filter housing and pipes. Touching these can cause severe burns. Therefore, let the plant cool down.
- Please make sure that the housing is in a pressureless state before opening it. If necessary, the housing can be depressurized.
- Steam will escape during the steaming procedure. Avoid any contact.

5. Scientific and Laboratory Services

Pall provides a full laboratory and field technical service to assist in the application and evaluation of Pall filter products. If you have technical questions please do not hesitate to use this customer service, available through your local sales office.

Figure 1: Recommended Filter Installation for in situ Steam Sterilization of Pall SUPRAdisc Modules

Bad Kreuznach – Germany

+49 671 8822 0 telephone

+49 671 8822 200 fax

Portsmouth – United Kingdom

+44 23 92 30 22 69 telephone

+44 23 92 30 25 09 fax

New York – USA

+1 516 484 5400 telephone

+1 516 625 3610 fax

Melbourne – Australia

+61 395 8481 00 telephone

+61 395 8466 47 fax

Paris – France

+33 1 30 6138 00 telephone

+33 1 30 6157 08 fax

Milano – Italy

+39 02 47 79 61 telephone

+39 02 41 22 985 fax

Filtration. Separation. Solution.sm

Visit us on the Web at www.pall.com

Pall Corporation has offices and plants throughout the world in locations including: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, United Kingdom, United States, and Venezuela. Distributors are located in all major industrial areas of the world.

© Copyright 2018, Pall Corporation. Pall,  and SUPRAdisc are trademarks of Pall Corporation.

® Indicates a Pall trademark registered in the USA. *Filtration. Separation. Solution.sm* is a service mark of Pall Corporation.

Because of developments in technology these data or procedures may be subject to change. Consequently we advise users to review their continuing validity annually.

Reorder Code. PPC-G125 dt./engl./franz./ital./span./jap./russ.